

## Beta-galactoside alpha-2,6-sialyltransferase 1

Catalog No: #AP88557



Package Size: #AP88557-1 50ug #AP88557-2 100ug #AP88557-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-galactoside alpha-2,6-sialyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Alpha 2,6-ST 1,CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,6-sialyltransferase 1,ST6Gal I,ST6GalI,Sialyltransferase 1,SIAT1
Accession No.	Uniprot ID:Q92182
Uniprot	Q92182
GeneID	396169;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ST6GAL1

Note: This product is for in vitro research use only