

NAD-dependent protein deacetylase sirtuin-7

Catalog No: #AP88558



Package Size: #AP88558-1 50ug #AP88558-2 100ug #AP88558-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent protein deacetylase sirtuin-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SIR2-like protein 7,SIR2L7
Accession No.	Uniprot ID:Q9NRC8
Uniprot	Q9NRC8
GeneID	51547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SIRT7

Note: This product is for in vitro research use only