

Serine racemase

Catalog No: #AP88561



Package Size: #AP88561-1 50ug #AP88561-2 100ug #AP88561-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine racemase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	D-serine ammonia-lyase,D-serine dehydratase,L-serine ammonia-lyase,L-serine dehydratase
Accession No.	Uniprot ID:Q9GZT4
Uniprot	Q9GZT4
GeneID	63826;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRR

Note: This product is for in vitro research use only