

## Heterogeneous nuclear ribonucleoproteins A2/B1

Catalog No: #AP88563



Package Size: #AP88563-1 50ug #AP88563-2 100ug #AP88563-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Heterogeneous nuclear ribonucleoproteins A2/B1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	hnRNP A2/B1,Hnrnp,Hnrpa2b1
Accession No.	Uniprot ID:A7VJC2
Uniprot	A7VJC2
GeneID	362361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hnrnpa2b1

Note: This product is for in vitro research use only