

NEDD8-activating enzyme E1 catalytic subunit

Catalog No: #AP88568



Package Size: #AP88568-1 50ug #AP88568-2 100ug #AP88568-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NEDD8-activating enzyme E1 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NEDD8-activating enzyme E1C, Ubiquitin-activating enzyme E1C, Ubiquitin-like modifier-activating enzyme 3, Ubiquitin-activating enzyme 3, UBE1C
Accession No.	Uniprot ID:Q8TBC4
Uniprot	Q8TBC4
GeneID	9039;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: UBA3

Note: This product is for in vitro research use only