

Urocortin-3

Catalog No: #AP88570



Package Size: #AP88570-1 50ug #AP88570-2 100ug #AP88570-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Urocortin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Stresscopin,Urocortin III,Ucn III,SPC
Accession No.	Uniprot ID:Q969E3
Uniprot	Q969E3
GeneID	114131;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UCN3

Note: This product is for in vitro research use only