

Serine/threonine-protein kinase TBK1

Catalog No: #AP88575



Package Size: #AP88575-1 50ug #AP88575-2 100ug #AP88575-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase TBK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NF-kappa-B-activating kinase,T2K,TANK-binding kinase 1,NAK
Accession No.	Uniprot ID:Q9UHD2
Uniprot	Q9UHD2
GeneID	29110;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TBK1

Note: This product is for in vitro research use only