

Toll-like receptor 9

Catalog No: #AP88578



Package Size: #AP88578-1 50ug #AP88578-2 100ug #AP88578-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toll-like receptor 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CD289
Accession No.	Uniprot ID:Q5I2M8
Uniprot	Q5I2M8
GeneID	403502;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TLR9

Note: This product is for in vitro research use only