

Thyrotropin receptor

Catalog No: #AP88581



Package Size: #AP88581-1 50ug #AP88581-2 100ug #AP88581-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thyrotropin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Thyroid-stimulating hormone receptor,TSH-R
Accession No.	Uniprot ID:P47750
Uniprot	P47750
GeneID	22095;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tshr

Note: This product is for in vitro research use only