

# Trefoil factor 1

Catalog No: #AP88582



Package Size: #AP88582-1 50ug #AP88582-2 100ug #AP88582-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Trefoil factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Breast cancer estrogen-inducible protein,PNR-2,Polypeptide P1.A,hP1.A,Protein pS2,BCEI,PS2
Accession No.	Uniprot ID:P04155
Uniprot	P04155
GeneID	7031;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TFF1

Note: This product is for in vitro research use only