

Trefoil factor 1

Catalog No: #AP88583



Package Size: #AP88583-1 50ug #AP88583-2 100ug #AP88583-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trefoil factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein pS2,Bcei,Ps2
Accession No.	Uniprot ID:Q08423
Uniprot	Q08423
GeneID	21784;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tff1

Note: This product is for in vitro research use only