

Mitochondrial brown fat uncoupling protein 1

Catalog No: #AP88586



Package Size: #AP88586-1 50ug #AP88586-2 100ug #AP88586-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitochondrial brown fat uncoupling protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	UCP 1,Solute carrier family 25 member 7,Thermogenin,Slc25a7,Ucp
Accession No.	Uniprot ID:P12242
Uniprot	P12242
GeneID	22227;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ucp1

Note: This product is for in vitro research use only