

Uveal autoantigen with coiled-coil domains and ankyrin repeats

Catalog No: #AP88588

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88588-1 50ug #AP88588-2 100ug #AP88588-3 1mg

Description

Product Name	Uveal autoantigen with coiled-coil domains and ankyrin repeats
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KIAA1561
Accession No.	Uniprot ID:Q9BZF9
Uniprot	Q9BZF9
GeneID	55075;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UACA

Note: This product is for in vitro research use only