

Succinate-semialdehyde dehydrogenase, mitochondrial

Catalog No: #AP88589

Package Size: #AP88589-1 50ug #AP88589-2 100ug #AP88589-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Succinate-semialdehyde dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Aldehyde dehydrogenase family 5 member A1,NAD(+)-dependent succinic semialdehyde dehydrogenase,SSADH
Accession No.	Uniprot ID:P51649
Uniprot	P51649
GeneID	7915;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALDH5A1

Note: This product is for in vitro research use only