

Estrogen sulfotransferase, testis isoform

Catalog No: #AP88590



Package Size: #AP88590-1 50ug #AP88590-2 100ug #AP88590-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Estrogen sulfotransferase, testis isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Sulfotransferase 1E1,ST1E1,Sulfotransferase, estrogen-preferring,Ste
Accession No.	Uniprot ID:P49891
Uniprot	P49891
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sult1e1

Note: This product is for in vitro research use only