

Syntaxin-1A

Catalog No: #AP88592



Package Size: #AP88592-1 50ug #AP88592-2 100ug #AP88592-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syntaxin-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Neuron-specific antigen HPC-1, Synaptotagmin-associated 35 kDa protein, P35A, Sap
Accession No.	Uniprot ID:P32851
Uniprot	P32851
GeneID	116470;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Stx1a

Note: This product is for in vitro research use only