

Prostaglandin E synthase 3

Catalog No: #AP88593



Package Size: #AP88593-1 50ug #AP88593-2 100ug #AP88593-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin E synthase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytosolic prostaglandin E2 synthase,cPGES,Hsp90 co-chaperone,Progesterone receptor complex p23,Telomerase-binding protein p23,P23,TEBP
Accession No.	Uniprot ID:Q15185
Uniprot	Q15185
GeneID	10728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTGES3

Note: This product is for in vitro research use only