

Tumor necrosis factor ligand superfamily member 14

Catalog No: #AP88594



Package Size: #AP88594-1 50ug #AP88594-2 100ug #AP88594-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Herpes virus entry mediator ligand,HVEM-L,CD258,HVEML,LIGHT,UNQ391/PRO726
Accession No.	Uniprot ID:O43557
Uniprot	O43557
GeneID	8740;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFSF14

Note: This product is for in vitro research use only