

Tumor necrosis factor ligand superfamily member 8

Catalog No: #AP88595



Package Size: #AP88595-1 50ug #AP88595-2 100ug #AP88595-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD30 ligand,CD30-L,CD153,CD30L,CD30LG
Accession No.	Uniprot ID:P32971
Uniprot	P32971
GeneID	944;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFSF8

Note: This product is for in vitro research use only