

THAP domain-containing protein 4

Catalog No: #AP88599



Package Size: #AP88599-1 50ug #AP88599-2 100ug #AP88599-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	THAP domain-containing protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CGI-36,PP238
Accession No.	Uniprot ID:Q8WY91
Uniprot	Q8WY91
GeneID	51078;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: THAP4

Note: This product is for in vitro research use only