

Alpha-globin transcription factor CP2

Catalog No: #AP88603



Package Size: #AP88603-1 50ug #AP88603-2 100ug #AP88603-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-globin transcription factor CP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SAA3 enhancer factor, Transcription factor LSF, LSF, SEF
Accession No.	Uniprot ID:Q12800
Uniprot	Q12800
GeneID	7024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TF_{CP2}

Note: This product is for in vitro research use only