

Tumor necrosis factor receptor superfamily member 12A

Catalog No: #AP88604

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88604-1 50ug #AP88604-2 100ug #AP88604-3 1mg

Description

Product Name	Tumor necrosis factor receptor superfamily member 12A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fibroblast growth factor-inducible immediate-early response protein 14,FGF-inducible 14,Tweak-receptor,TweakR,CD266,FN14
Accession No.	Uniprot ID:Q9NP84
Uniprot	Q9NP84
GeneID	51330;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFRSF12A

Note: This product is for in vitro research use only