

General transcription factor IIH subunit 2

Catalog No: #AP88608



Package Size: #AP88608-1 50ug #AP88608-2 100ug #AP88608-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	General transcription factor IIH subunit 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	General transcription factor IIH polypeptide 2
Accession No.	Uniprot ID:Q2TBV5
Uniprot	Q2TBV5
GeneID	511829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GTF2H2

Note: This product is for in vitro research use only