

## Thyroid hormone receptor beta

Catalog No: #AP88611



Package Size: #AP88611-1 50ug #AP88611-2 100ug #AP88611-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thyroid hormone receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nuclear receptor subfamily 1 group A member 2,c-erbA-2,c-erbA-beta,Erba2,Nr1a2
Accession No.	Uniprot ID:P37242
Uniprot	P37242
GeneID	21834;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Thrb

Note: This product is for in vitro research use only