

Urocortin-2

Catalog No: #AP88612



Package Size: #AP88612-1 50ug #AP88612-2 100ug #AP88612-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Urocortin-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Stresscopin-related peptide,Urocortin II,Ucn II,Urocortin-related peptide,SRP,URP |
| Accession No. | Uniprot ID:Q96RP3 |
| Uniprot | Q96RP3 |
| GeneID | 90226; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UCN2

Note: This product is for in vitro research use only