

## Sortilin-related receptor

Catalog No: #AP88614



Package Size: #AP88614-1 50ug #AP88614-2 100ug #AP88614-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sortilin-related receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gp250,Low-density lipoprotein receptor relative with 11 ligand-binding repeats,LDLR relative with 11 ligand-binding repeats,SorLA-1,Sorting protein-related receptor containing LDLR class A repeats,mSorLA
Accession No.	Uniprot ID:O88307
Uniprot	O88307
GeneID	20660;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sort1

Note: This product is for in vitro research use only