

## Serine/threonine-protein kinase TAO1

Catalog No: #AP88617



Package Size: #AP88617-1 50ug #AP88617-2 100ug #AP88617-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase TAO1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kinase from chicken homolog B,hKFC-B,MARK Kinase,MARKK,Prostate-derived sterile 20-like kinase 2,PSK-2,Thousand and one amino acid protein kinase 1,TAOK1,KIAA1361,MAP3K16,MARKK
Accession No.	Uniprot ID:Q7L7X3
Uniprot	Q7L7X3
GeneID	57551;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TAOK1

Note: This product is for in vitro research use only