

Transforming growth factor beta-1

Catalog No: #AP88622



Package Size: #AP88622-1 50ug #AP88622-2 100ug #AP88622-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TGF-beta-1
Accession No.	Uniprot ID:P18341
Uniprot	P18341
GeneID	282089;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TGFB1

Note: This product is for in vitro research use only