

Succinate-semialdehyde dehydrogenase, mitochondrial

Catalog No: #AP88626

Package Size: #AP88626-1 50ug #AP88626-2 100ug #AP88626-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Succinate-semialdehyde dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aldehyde dehydrogenase family 5 member A1,NAD(+)-dependent succinic semialdehyde dehydrogenase
Accession No.	Uniprot ID:Q8BWF0
Uniprot	Q8BWF0
GeneID	214579;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh5a1

Note: This product is for in vitro research use only