

Bile salt sulfotransferase

Catalog No: #AP88627



Package Size: #AP88627-1 50ug #AP88627-2 100ug #AP88627-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile salt sulfotransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hydroxysteroid sulfotransferase,ST,ST-20,Sulfotransferase 2A1,ST2A1
Accession No.	Uniprot ID:P15709
Uniprot	P15709
GeneID	24912;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:St2a1

Note: This product is for in vitro research use only