

Signal recognition particle subunit SRP68

Catalog No: #AP88629



Package Size: #AP88629-1 50ug #AP88629-2 100ug #AP88629-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Signal recognition particle subunit SRP68
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SRP68,Signal recognition particle 68 kDa protein
Accession No.	Uniprot ID:Q8BMA6
Uniprot	Q8BMA6
GeneID	217337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Srp68

Note: This product is for in vitro research use only