

Tumor necrosis factor ligand superfamily member 11

Catalog No: #AP88630



Package Size: #AP88630-1 50ug #AP88630-2 100ug #AP88630-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Osteoclast differentiation factor,ODF,Osteoprotegerin ligand,OPGL,Receptor activator of nuclear factor kappa-B ligand,RANKL,TNF-related activation-induced cytokine,TRANCE,CD254,Opgl,Rankl,Trance
Accession No.	Uniprot ID:O35235
Uniprot	O35235
GeneID	21943;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfsf11

Note: This product is for in vitro research use only