

# Tumor necrosis factor receptor superfamily member 10A

Catalog No: #AP88632

Package Size: #AP88632-1 50ug #AP88632-2 100ug #AP88632-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 10A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Death receptor 4, TNF-related apoptosis-inducing ligand receptor 1, TRAIL receptor 1, CD261, APO2, DR4, TRAILR1
Accession No.	Uniprot ID:O00220
Uniprot	O00220
GeneID	8797;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TNFRSF10A

Note: This product is for in vitro research use only