

Tumor necrosis factor receptor superfamily member 11B

Catalog No: #AP88633

Package Size: #AP88633-1 50ug #AP88633-2 100ug #AP88633-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 11B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Osteoprotegerin,Opg
Accession No.	Uniprot ID:O08727
Uniprot	O08727
GeneID	25341;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfrsf11b

Note: This product is for in vitro research use only