

Somatostatin receptor type 5

Catalog No: #AP88638



Package Size: #AP88638-1 50ug #AP88638-2 100ug #AP88638-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Somatostatin receptor type 5 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | SS-5-R |
| Accession No. | Uniprot ID:P35346 |
| Uniprot | P35346 |
| GeneID | 6755; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SSTR5

Note: This product is for in vitro research use only