

## Tumor necrosis factor receptor superfamily member 9

Catalog No: #AP88641



Package Size: #AP88641-1 50ug #AP88641-2 100ug #AP88641-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Tumor necrosis factor receptor superfamily member 9                  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Ms   |
| Other Names        | 4-1BB ligand receptor, T-cell antigen 4-1BB, CD137, Cd137, Ila, Ly63 |
| Accession No.      | Uniprot ID:P20334  |
| Uniprot            | P20334   |
| GeneID             | 21942;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)               |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Tnfrsf9

Note: This product is for in vitro research use only