

Bile salt sulfotransferase

Catalog No: #AP88644



Package Size: #AP88644-1 50ug #AP88644-2 100ug #AP88644-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile salt sulfotransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dehydroepiandrosterone sulfotransferase,DHEA-ST,Hydroxysteroid Sulfotransferase,HST,ST2,ST2A3,Sulfotransferase 2A1,ST2A1,HST,STD
Accession No.	Uniprot ID:Q06520
Uniprot	Q06520
GeneID	6822;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SULT2A1

Note: This product is for in vitro research use only