

# Toll-like receptor 4

Catalog No: #AP88645



Package Size: #AP88645-1 50ug #AP88645-2 100ug #AP88645-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Toll-like receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD284
Accession No.	Uniprot ID:Q9GL65
Uniprot	Q9GL65
GeneID	281536;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TLR4

Note: This product is for in vitro research use only