

Uridine phosphorylase 1

Catalog No: #AP88652



Package Size: #AP88652-1 50ug #AP88652-2 100ug #AP88652-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Uridine phosphorylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	UPase 1,UP
Accession No.	Uniprot ID:Q16831
Uniprot	Q16831
GeneID	7378;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UPP1

Note: This product is for in vitro research use only