

Ubiquitin-conjugating enzyme E2 B

Catalog No: #AP88656



Package Size: #AP88656-1 50ug #AP88656-2 100ug #AP88656-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-conjugating enzyme E2 B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	E2 ubiquitin-conjugating enzyme B,RAD6 homolog B,HR6B,Ubiquitin carrier protein B,Ubiquitin-conjugating enzyme E2-17 kDa,Ubiquitin-protein ligase B,RAD6B
Accession No.	Uniprot ID:P63146
Uniprot	P63146
GeneID	7320;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UBE2B

Note: This product is for in vitro research use only