

## Tenomodulin

Catalog No: #AP88662



Package Size: #AP88662-1 50ug #AP88662-2 100ug #AP88662-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tenomodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TeM, Chondromodulin-1-like protein, ChM1L, Chondromodulin-I-like protein, Myodulin, Tendin, Chm1l
Accession No.	Uniprot ID:Q9EP64
Uniprot	Q9EP64
GeneID	64103;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Tnmd

Note: This product is for in vitro research use only