

TGF-beta receptor type-2

Catalog No: #AP88664



Package Size: #AP88664-1 50ug #AP88664-2 100ug #AP88664-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | TGF-beta receptor type-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | TGFR-2,TGF-beta type II receptor,Transforming growth factor-beta receptor type II,TGF-beta receptor type II |
| Accession No. | Uniprot ID:P38551 |
| Uniprot | P38551 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TGFBR2

Note: This product is for in vitro research use only