

Stathmin

Catalog No: #AP88665



Package Size: #AP88665-1 50ug #AP88665-2 100ug #AP88665-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stathmin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Leukemia-associated phosphoprotein p18, Metablastin, Oncoprotein 18, Op18, Phosphoprotein p19, pp19, Pr22 protein, ProsoLin, pp17, Lap18
Accession No.	Uniprot ID:P13668
Uniprot	P13668
GeneID	29332;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Stmn1

Note: This product is for in vitro research use only