

Steryl-sulfatase

Catalog No: #AP88666



Package Size: #AP88666-1 50ug #AP88666-2 100ug #AP88666-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Steryl-sulfatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Arylsulfatase C,ASC,Steroid sulfatase,Steryl-sulfate sulfohydrolase
Accession No.	Uniprot ID:P15589
Uniprot	P15589
GeneID	24800;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sts

Note: This product is for in vitro research use only