

Transmembrane protein 158

Catalog No: #AP88668



Package Size: #AP88668-1 50ug #AP88668-2 100ug #AP88668-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transmembrane protein 158
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Ras-induced senescence protein 1, RIS1
Accession No.	Uniprot ID:A2VDX9
Uniprot	A2VDX9
GeneID	788085;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TMEM158

Note: This product is for in vitro research use only