

Tumor necrosis factor receptor superfamily member 12A

Catalog No: #AP88669

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88669-1 50ug #AP88669-2 100ug #AP88669-3 1mg

Description

Product Name	Tumor necrosis factor receptor superfamily member 12A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fibroblast growth factor-inducible immediate-early response protein 14,FGF-inducible 14,Fibroblast growth factor-regulated protein 2,Tweak-receptor,TweakR,CD266,Fgfrp2,Fn14
Accession No.	Uniprot ID:Q9CR75
Uniprot	Q9CR75
GeneID	27279;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfrsf12a

Note: This product is for in vitro research use only