

Pro-thyrotropin-releasing hormone

Catalog No: #AP88671



Package Size: #AP88671-1 50ug #AP88671-2 100ug #AP88671-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pro-thyrotropin-releasing hormone
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pro-TRH, Prothyroliberin
Accession No.	Uniprot ID:P20396
Uniprot	P20396
GeneID	7200;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TRH

Note: This product is for in vitro research use only