

TATA box-binding protein-like protein 1

Catalog No: #AP88672



Package Size: #AP88672-1 50ug #AP88672-2 100ug #AP88672-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TATA box-binding protein-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TBP-like protein 1,21 kDa TBP-like protein,Second TBP of unique DNA protein,STUD,TATA box-binding protein-related factor 2,TBP-related factor 2,TBP-like factor,TBP-related protein,TLF,TLP,TLP21,TRF2,TRP
Accession No.	Uniprot ID:P62380
Uniprot	P62380
GeneID	9519;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TBPL1

Note: This product is for in vitro research use only