

Transmembrane protein 158

Catalog No: #AP88678



Package Size: #AP88678-1 50ug #AP88678-2 100ug #AP88678-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transmembrane protein 158
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	40 kDa BINP-binding protein,p40BBP,Ras-induced senescence protein 1,Mbbp,Ris1
Accession No.	Uniprot ID:Q6F5E0
Uniprot	Q6F5E0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tmem158

Note: This product is for in vitro research use only