

Serine/arginine-rich splicing factor 1

Catalog No: #AP88684



Package Size: #AP88684-1 50ug #AP88684-2 100ug #AP88684-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/arginine-rich splicing factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Splicing factor, arginine/serine-rich 1,SFRS1,RCJMB04_1I5
Accession No.	Uniprot ID:Q5ZML3
Uniprot	Q5ZML3
GeneID	772264;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRSF1

Note: This product is for in vitro research use only